

LIGHT & BRIGHT™ TREATMENT

Ten things you should know now.

A recent industry survey underscores that markets for IPL and fractional resurfacing continue to grow.¹ Taking these two business-expanding technologies to the next level, Candela's Light & Bright treatment protocol unites complementary IPL and non-ablative fractional treatments to deliver high rates of improvement and even higher patient satisfaction.² The new Light & Bright treatment, available on the Nordlys™ multi-application platform, uses patented Ellipse IPL™ and the Frax 1550™ handpieces to address both vascular and pigmented lesions as well as textural irregularities via skin resurfacing.

10 THINGS YOU SHOULD KNOW ABOUT LIGHT & BRIGHT TREATMENT

1. Offers complementary Ellipse IPL™ and Frax 1550™ treatments on the Nordlys platform^{3,4}
2. Treats diffuse redness, pigmentation, and textural irregularities on the face and/or décolleté via skin resurfacing^{3,4}
3. Demonstrates impressive results in just three treatment sessions at one-month intervals^{3,4}
4. Ellipse IPL technology uses targeted (narrowband), controlled (sub-millisecond pulses), and filtered light (dual filtration technology)³⁻⁵
5. Ellipse IPL technology uses less than half the fluence, no active cooling, and requires fewer treatments compared to broadband IPL³⁻⁵
6. The Frax 1550 handpiece delivers high ROI treatments with a smart approach to minimizing consumables³
7. Treatment results in high improve-

Light & Bright™



PR 530,
VL 555
or PL 400



Frax 1550



Frax 1550 or
Ellipse IPL (PR 530,
VL 555 or PL 400)
depending on skin type

Facial and chest areas treated by the Ellipse IPL and Frax 1550 Handpieces


8. Treatment achieves high patient satisfaction—95% of treated patients report they are satisfied or extremely satisfied one month
9. Candela supports practices offering the Light & Bright treatment with practice marketing materials, digital ads, and geo-targeted social media boosting

ment rates—patients experience 80% improvement in telangiectasias and diffuse redness²

after treatment completion with improvement in skin texture* and in signs of photoaging^{2,**}

10. Your patients will love it! Treatment leaves skin lighter, brighter, and smoother²

95% SATISFIED AT 1-MONTH FOLLOW-UP VISIT (N=20)



CLINICAL STUDY RESULTS

A clinical study of 22 subjects found high patient satisfaction following a series of three Ellipse IPL and Frax 1550 treatments to the face and/or chest. Subject satisfaction continued to improve during the treatment series, particularly for those patients who were extremely satisfied at the one-month follow-up visit. Investigator assessments showed improvement in pigmentary and vascular dyschromia after the first treatment and continued to improve with successive treatments. Clinical photography demonstrated visible improvement in photoaged skin. The versatility of having access to a selection of Ellipse IPL handpieces with dual filtering and a fractional resurfacing handpiece in one multi-application, multi-technology platform allows for the treatment of various clinical indications in one treatment session.^{3,5}

A FOUNDATIONAL PRACTICE DEVICE: NORDLYS SYSTEM

Meet your patients' increasing demands for aesthetic and dermatological procedures with a single multi-appli-

PATIENTS LOVE THE LIGHT & BRIGHT TREATMENT



Photo credit: Photos courtesy of Konika Patel Schallen, MD

Diffuse redness and individual telangiectasia; textural irregularity.* Baseline (left) and 12 weeks post third treatment (right). Photos are unretouched. Individual results may vary.

"The Light and Bright treatment decreased my redness drastically, in very little time. I feel so much more confident and comfortable in my own skin. It has truly changed my life."

— AMY, LOWELL, MA



cation platform for your practice—the Nordlys System.^{3,4}

- 8 Ellipse IPL applicators
- Expandable fractional capabilities with Frax 1550 and Frax 1940⁷ applicators
- Nd:YAG 1064 applicator
- Extensive treatment options with 21 FDA cleared indications ranging

from hair removal to vascular and pigmentary concerns

- Quiet⁶ and ergonomic operation
- Secure, Wi-Fi-enabled remote assistance

Learn more about the Light & Bright treatment and the versatile Nordlys system platform at Candelamedical.com/Nordlys10Things ■

* Improvement of textural irregularities via skin resurfacing

** Via reduction in pigmentation, vascularity and improvement in textural irregularities via skin resurfacing.

1. The American Society for Aesthetic Plastic Surgery 2019 Cosmetic (Aesthetic) Surgery National Bank Statistics. <https://www.surgery.org/sites/default/files/ASAPS-Stats2018.pdf>. Accessed March 18, 2020.
 2. Light & Bright Treatment. Data on file, Candela 2020.
 3. Ellipse Nordlys 510(k) clearance (K161162), September 2016
 4. Nordlys system instructions for use, 2018. Candela, data on file.
 5. Bjerring P, Christiansen K, Troilius A, Dierickx C. Facial photo rejuvenation using two different intense pulsed light (IPL) wavelength bands. *Lasers Surg Med.* 34:120-126. doi:10.1002/lsm.20000
 6. Data on file, Candela, 2020.
 7. Ellipse Frax Pro and Nordlys 510(k) clearance (K192951), March 2020.

© 2020 Candela Corporation. This material contains registered and unregistered trademarks, trade-names, service marks and brand names of Candela Corporation and its affiliates. All other trademarks are the property of their respective owners. All rights reserved.